Bluetooth Connection and Data Syncing (1) Enable your phone's Bluetooth, open the Sleepus app, and turn on the device. Follow the instructions on the app to bind pair and connect the device. (2) Open the app and, tap press the device button on the device. The Bluetooth light will illuminate and turn on, and the data will be synced to the app.

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This product is a physiotherapy massager which is designed to massage and vibrate and stimulate the respiratory vestibular nerve in the ear muscles to keep the respiratory tract airway open when during sleeping.

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Wearing Instructions Rotate and use the product Aas shown below. Adjust the product to a stable and comfortable position.

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
Power On/Off: Press and hold the power switch for 3 seconds to turn on the device. The power indicator will automatically turn off after 10 seconds. Press the power switch again to turn on the power indicator on again. To turn the device off, simply pPress and hold the power switch for 3 seconds to turn off the device.

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Adjusting the intensity of intervention: If the intensity of stimulation is too high, you can log in to the app and enter the Settings interface to adjust the intensity of intervention.

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